Combi

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
# Index

- Precautions ..... 1
- Name of Components ..... 2
- Cleaning and Disinfection ..... 3
- Preparation before Breast Milk Expression ..... 5
- Double Electric Breast Pump Installation ..... 7
- Double Electric Breast Pump Motor Descriptions ..... 9
- Power Connection Method ..... 10
- Double Electric Breast Pump Operation ..... 11
- Convert Electric Breast Pump to Manual Type ..... 12
- Storing and Preparing Breastmilk ..... 13
- Combi Breast Feeding Support Series..... 14
- Troubleshooting ..... 15
- Product Specifications ..... 16

## Precautions

- Read this user manual in detail before use. Check there is no cracks or damages on any of the parts.
- Clean your hands with detergent before use.
- Clean all accessories (excluding the main body and silicone tube) with Combi Detergent for Feeding Bottles. After that, boil or steam for sterilization. If sterilized by boiling water, the recommended duration is 2 ~ 3 minutes, but never exceeds 3 minutes. Stir the boiling water after putting in the accessories in order to completely sterilize them in the flowing water.
- Keep the plastic accessories away from direct sunlight after clean in order to prevent aging components.
- Atomization after sterilization is normal phenomenon for silicone accessories due to its water absorption characteristics.
- Long-term use, scrub, sterilization and drying of the plastic accessories may cause deformation or deterioration which is unavoidable. Replace the accessory when it happens.
- Do not tilt the breast milk storage bottle or store over the maximum capacity line.
- In case the milk storage bottle is tilted, switch off the pump and then correct the position.
- Consult and use the pump as advised by doctor or nurse for users with hollow nipple, mastitis, flat nipple or inadequate breast milk produce after giving birth.
- Do not dismantle, repair or alter the main body on your own.

### Warnings

- Do not use the pump when there is crack or damage on the power cord or when the plug is loose.
- Do not pull or lick the power cord. Keep the power cord away from children when it is not used.
- Remove the power cord from the socket if the main body gets wet immediately and dry the main body.
- Remove the power cord when it is not used and keep the pump and all the accessories away from children.
- Ensure there is nothing attached on the power cord and on the connector before connecting the power cord to avoid short circuit.
- Rub away dust on the power cord regularly in order to avoid short circuit.
- Hold the plug end and remove the power. Do not pull the power cord upon power removal.
- Do not touch the power cord, plug or socket when your hands are not dry in order to avoid the risk of electric shock.
- Combi electrical breast pump accessories are exclusive for this pump. Do not use accessories of the other brands in order to avoid poor pumping effects and damages of the main body.
- Do not install glass milk bottle onto this pump in order to avoid accidents.
- Attached power cord is exclusive for this pump and not applicable to other electronic products.
- People with insufficient fitness, perception, awareness, no experience and knowledge (including children) should not use the pump, unless under guardian supervision.

## Name of Components

### Electric Breast Pump

### Accessories

## Clean and Disinfection

Cleaning before the first use and after every usage

**1** Clean all accessories (except the main body, adaptor, silicone tube, T-shaped tube, and USB cable).

**2** Rinse all accessories in warm water.

**3** Completely clean all accessories with bottle brush, nipple brush and Combi Detergent for Feeding Bottles.

**4** Carefully rinse all accessories with warm water.

**5-1** Sterilize in boiling water.

**5-2** Or sterilize in steam sterilizer.

- 6 After sterilization, place them onto a dry towel or put them into the storage box directly.

### Notes

- Clean all accessories with warm water and detergent to avoid milk scale and bacteria, expect main body, silicone tube, adaptor, T-shaped tube and USB cable.
- Do not use hard bristles or abrasive powder for cleaning to avoid damaging the products. Clean with soft sponge.
- Replace the valve when it is damaged or its adaptation is poor.
- Rinse the milk attached on the valve in warm water before cleaning the pump and the valve, so as to remove the valve easier without fat in the milk. Excessive force will damage the valve.
- Sterilize accessories by boiling water or sterilizer after cleaning. If sterilized by boiling water, the recommended duration is 2 to 3 minutes, but never exceeds 3 minutes. Stir the boiling water after putting in the parts in order to completely sterilize them in the flowing water and prevent deformation.
- Air-dry accessories in a cool dry location after sterilization. Ensure there is no any water remained before setting up the pump.
- Do not immerse or rinse the main body and the adapter in water. Rub with a clean wet towel if dirty.
- Avoid excessive force when unloading accessories for cleaning or may cause damages to the accessories.
- Keep all accessories properly to avoid missing.
- Do not wash this product with automatic dishwasher.
- Do not clean this product with bleach.

## Preparation before Breast Milk Expression

### Massage before expressing

Step 1. Having a happy mood which is very important to produce pleasant breast milk for our babies.

Step 2. Clean your hands and complete the disinfection steps of accessories in order to protect the baby.

Step 3. Cleaning nipple. However, in general, your nipple does not require any special cleaning unless you are sweating. You should avoid over cleaning nipple with detergent or soap which may result dry and cracked nipple.

Step 4. Breast massage. In general, no preparation is needed before pumping if milk runs out smoothly and there is no lump in the breast. However, if there is lump, you may massage your breast timely. First, press 2 to 3 cm from the nipple to stimulate milk ejection. Never massage the breast firmly or strongly to avoid injuring breast.


#### Massage method

**Circular Massage:**  
Massage the breast back and forth and up and down with both hands.

**Spiral Massage:**  
Fix the breast with one hand and massage from base to nipple with another hand orderly.

**Extrusion Massage:**  
Hold both sides of the breast with both hands, press from base to nipple.

Step 5. Pumping on alternating sides until milk flow slows down, which can stimulate the milk expression to increase milk flow speed.

## Warm Paste? Cold Paste? What are the differences?

Warm paste can stimulate milk secretion which is suitable for mother with soft breast after production and not yet with breast tenderness, or with small milk amount. Be careful the temperature should not exceed 40°C to avoid skin burn. Do not cover nipple. Warm paste should limit to 10 to 15 minutes but not too long. You may also gently lean your body in bath. The steam during bath can also result in the same effect.

However, warm paste is not suitable for stuffed breast, swelling breast, or pain breast after breast milk expression. Cold Paste should be used instead in order to reduce breast surface temperature which reduces pain, edema and inflammation reactions.

Use water bag, wet towels or washed cabbage leaves. Do not use ice directly to avoid frostbite due to too low temperature.

Do not cold paste areola and nipple, otherwise it will suppress milk discharge.

Discharge breast milk correctly is the most important regardless warm or cold paste.

Consult doctor immediately if your breast skin is pain, swelling, fewer milk or you feel cold and get fever, but still continue emptying your breast.

### Warm Paste

### Cold Paste

## Double Electric Breast Pump Installation

### Setting up Double Electric Breast Pump

- 1 Insert the round bottom of valve into the circular groove of the bottom of shield breast.

! To facilitate disassemble, 1cm of insertion is suggested.

- 2 Connect and lock the breast shield with the feeding bottle to avoid leaking.

- 3 Insert the silicone diaphragm into the circular groove on top of the breast shield.

! Ensure diaphragm edge has securely inserted into the circular groove.

- 4 Insert the silicone cover into the funnel of breast shield.

- 5 Cover silicone diaphragm with top cover.

! Remove the top cover from bottom to top during disassembly.

### 6-1 Using Single Breast Pump

**Method 1 :** To start using, connect the silicone tube with the top cover and air valve of the main body. Connect power adaptor with main body and power socket (or using USB power battery).

**Method 2 :** To start using, insert T-shaped tube into the air valve of the main body. Close one end of the T-shaped tube. Connect the other end with top cover and connect power adaptor with main body and power socket (or using USB power battery).

### 6-2 Using Double Breast Pump

To start using, insert T-shaped tube into the air valve of the main body and connect two silicone tubes with the T-shaped tube at one end and both top covers at the other. Then, connect the power adaptor with main body and power socket (or using USB power battery).

**Note:** Make sure the body is set on correct mode when using Single/Double function to avoid insufficient or excessive suction power.

**Double mode:** Please use two sets of breast shields.

**Single mode:** Please use one set of breast shield only.

(Make sure the other end of the T-shaped tube is closed when using.)

## Double Electric Breast Pump Motor Descriptions

Screen display / button	Description
Power Button	<p>Power on : After connecting with the power source, press the button and hold for 3sec. to start up the body. After the machine is on, press the button one time to start collecting milk.</p> <p>Pause : During the use, you can press the button at any time to pause the milk collection. After pausing for 2mins, the machine will shut down automatically to save the power. Press the power button again can restart milk collecting function.</p> <p>Power off : Press the button and hold for 3sec. during the use can power off the machine.</p>
Simulation / Expression Mode Switching Key	Simulation is the default mode once the power is on. Running with this default mode for two minutes the machine will automatically switch to milk expression mode. You can switch between these two modes during the use (The machine memorizes the last simulation / expression power level while reset is not necessary.)
Increase Power Level Single/Double Mode Switch	To increase power during use. Highest level is 10. Press and hold the button can switch the mode from single to double.
Decrease Power Level Single/Double Mode Switch	To decrease power during use. Lowest level is 1. Press and hold the button can switch the mode from double to single.
Simulation / Expression Mode	<ul style="list-style-type: none"> <li>💧 Simulation : Breast massage to stimulate lactation.</li> <li>💧 Expression : Milk expressing and deep pumping.</li> </ul>
Power Level	To display current suction power level.
Simulation / Expression Timer	To record using time after start and keep recording when mode switched.
Single/Double Mode Display	<ul style="list-style-type: none"> <li>🌸 Single-layer dynamic petal : Single breast pump mode.</li> <li>🌸 Double-layer dynamic petal : Double breast pump mode.</li> </ul>

## Power connection method

### 1 Using AC Power Adaptor:

Plug the L-connector end of the AC adaptor into the DC socket of the main body and insert the adaptor into the power socket.

### 2 Using USB Power Battery:

Plug the L-connector of USB cable into the main body and connect the USB plug with the USB power battery.  
(Power bank: 5V / 2.1A)

**Note:** Using AC Power Adaptor will best perform the product.

### Warnings

- Make sure the voltage and current are correct when using USB power battery. Low voltage will lead to insufficient power for the machine and over high voltage may cause damage to the machine.

Suggested Voltage: 5V~6V

Suggested Current: 1A~2.1A

\*Out of this range may fail to supply power for the product.

- To prevent foreign parts being accidentally sucked into the product and cause damage, make sure all parts of the product have been properly assembled before start using the device.
- To prevent potential risk of electric shock, do not connect the adaptor or USB cable to the power with wet hands.
- Notice the power supply of the local country and equip the product with suitable adaptor, if you use the product abroad.
- Use the adaptor and USB cable safely and keep them away from children.
- Normal operation temperature of this product ranges between 5°C~45°C.
- It is recommended not to use the product for more than 50 minutes and the interval between two-time uses shall be at least 40 minutes.

## Double Electric Breast Pump Operation

### Single Sided Suction

- Place the nipple at the center of the silicone cover. Fit the rim of the funnel tightly onto the breast. Grip the neck of the breast shield between the thumb and the index finger. Grab the bottle tightly with the rest of the fingers.
- Press the breast pump slightly against the breast at the start of the suction to ensure that the breast and the areola are fit tightly onto the funnel and achieve a vacuum effect required for suction.
- Adjust the level of the suction power according to personal comfort. Discomfort on the breast will occur if the power is too strong. If the power is too weak, the suction result will be affected.

### Double Sided Suction

- Hold a set of suction bottles on each hand and place the nipple at the center of each silicone cover. Fit the rim of the funnel tightly onto the breast. Grip the neck of the breast shield between the thumb and the index finger. Grab the bottle tightly with the rest of the fingers and start the suction.
- Adjust the level of the suction power according to personal comfort. Discomfort on the breast will occur if the power is too strong. If the power is too weak, the suction result will be affected.

### Precautions of using breast pump

- Read "Preparation before Breast Milk Expression" in this manual in order to enhance pumping effectiveness.
- Ensure your breast is tightly fit into the pump for vacuum pumping effect.
- "Bang" sound during use is air squeak sound, meaning that breast shield and breast is not closely attached.
- Select your most comfortable power level and speed during use, breast discomfort may be resulted if it is too strong, while pumping effect may be affected if it is too weak.
- Stop using the pump and consult professional doctors and nurses for discomfort during breast pumping.
- Do not store over the maximum capacity line of the breast milk storage bottle.
- The thread of the given breast shield is exclusively designed for Combi standard feeding bottle.
- Do not use feeding bottle of other brands to avoid milk leakage caused by poor adaptation.
- Do not use glass feeding bottle to avoid dropping the pump due to overweight after filling in breast milk.
- Always place the pump vertically to avoid leaking from breast shield or stoma.
- Replace a new white valve when it is damaged or its adaptation is poor since the adaptation of duckbill valve effects suction power.

## Convert Electric Breast Pump to Manual Type

*(Accessories of manual breast pump purchased separately)*

- 1 Assemble the silicone diaphragm and shaft.
- 2 Detach the tube, top cover and electric silicone diaphragm. Insert the manual silicone diaphragm into the circular groove on top of the breast shield.
- 3 Insert the U-shaped end of the handle into the shaft. Press down the handle to lock it with the breast shield. Place the top cover.

Ensure shaft and diaphragm installed securely.

Ensure diaphragm edge has securely inserted into the circular groove on the top.

Make sure a "click" sound is heard to assemble correctly.

### Notes for using the pump

- Please refer to the user manual of Combi Manual Breast Pump for detailed operation of the manual pump.
- Color of manual pump accessory (top cover, handle) is purple which is different from the electric breast bump but does not affect its functionality.

## Storing and Preparing Breastmilk

### *Freshly expressed breastmilk storage guidelines*

#### Breast Milk in Refrigeration

Storage Location	Environmental Temperature	Storage Time Limit	Container	Remarks
Room Temperature	below 25°C	6-8 hours	Milk Bottle	<ul style="list-style-type: none"> <li>• Use storage cap or cover.</li> <li>• Notice the temperature changes.</li> </ul>
Fridge	-15~4°C	24 hours	Milk bottle or storage bag	<ul style="list-style-type: none"> <li>• Use ice bag or ice brick.</li> <li>• Do not open the ice bucket too often.</li> </ul>
Freezer	0~4°C	5 days	Milk bottle or storage bag	<ul style="list-style-type: none"> <li>• Place in the inner part of the fridge.</li> <li>• Do not place in the outer part or near the door.</li> </ul>

#### Freezing Breast Milk

Storage Location	Environmental Temperature	Storage Time Limit	Container	Remarks
Freezer in single door fridge	-15°C	2 weeks	Milk Storage Bag	<ul style="list-style-type: none"> <li>• Using water-proof sticker to record date and volume of milk storage.</li> </ul>
Freezer in multiple doors fridge	-18°C	3 to 6 months	Milk Storage Bag	<ul style="list-style-type: none"> <li>• Do not store more than the maximum capacity stated on the storage bag to avoid extrusion split due to freezing expansion.</li> </ul>
Standalone Freezer	-20°C	6 to 12 months	Milk Storage Bag	<ul style="list-style-type: none"> <li>• Place the milk bottles or storage bags at the inner side of the refrigerator.</li> <li>• Feed on first-in-first-out principle.</li> </ul>

#### Thawing / Heating Breast Milk

State of Breast Milk	Thawing / Heating Method	Deposited time in room temperature after thaw	Breast Milk after thawing	Remarks
Fresh breast milk	Feed directly / Impermeable heating	X	Feed instantly	<ul style="list-style-type: none"> <li>• Do not re-freeze thawed breast milk.</li> <li>• Do not re-heat thawed breast milk.</li> </ul>
Taken out from fridge	Impermeable heating	X	Feed instantly	<ul style="list-style-type: none"> <li>• Do not use microwave oven to thaw or heat breastmilk.</li> </ul>
Taken out from freezer	Thaw in running cool water and then by impermeable heating	4 hours	Feed instantly	<ul style="list-style-type: none"> <li>• Impermeable heat in 60°C or below water.</li> <li>• Shake container softly while heating to mix the fat and milk well and ensure evenly heated.</li> </ul>
	Thaw in running cool water and then by impermeable heating	4 hours	To be consumed within 24 hours	

Source: 1.Taiwan National Health Council Year 2014 Breastfeeding Handbook.  
 2.The Academy of Breastfeeding Medicine (2004)Clinical Protocol Number #8 : Human Milk Storage Information for Home Use for Healthy Full Term Infants

## Combi Breast Feeding Support Series

Combi Breast Feeding Support

### Feeding

|| Combi Breast Feeding Bottles

### Maternity Underwear

|| Cupro underwear series  
|| Natural Colored Cotton underwear series

### Cleaning & Disinfection

|| High efficiency disinfection sterilizer and dryer.  
|| Milk bottles storage container

### Milk Collection Assistance

|| Combi Natural Rhyme Double Electric Breast Pump  
|| Combi Natural Rhyme Electric & Manual 2-in-1 Breast Pump  
|| Combi Natural Rhyme Electric Breast Pump  
|| Combi Natural Rhyme Manual Breast Pump

### Personal Care

|| 3D Breast Pad

## Troubleshooting

### Why is there no suction or insufficient suction power?

- The valve is not closely inserted into the breast shield or it is damaged.
- Silicone diaphragm is not closely inserted into the breast shield or it is damaged.
- Silicone tube is not closely inserted into the main body and the top cover or it is damaged.
- The top cover of the electric breast pump is not fitted tightly on to the top of breast shield.
- The power supply is not connected or the power from the portable power battery is not enough.
- The USB power supply is not reaching 5~6V.

### Why is milk overflowing or spilling out from the breast pump?

- Hold the breast pump vertically when in use.
- Please use the feeding bottle base to prevent the breast pump from overturning.
- The thread of the breast shield is exclusively designed for Combi standard feeding bottle. Do not use feeding bottle of other brands to avoid milk overflow caused by poor adaptation.

### Why is there breast milk residues remaining in the silicone diaphragm or silicone cover?

- Due to the amount of breast milk and strength of the suction power, breast milk will be ejected into the breast shield during pumping.
- Breast milk remained in silicone diaphragm or silicone cover is normal phenomenon but not caused by adaptation defect of the product.

### Why are silicone accessory parts become misty?

- Due to the water absorption property of silicone, it is normal for the silicone accessories to appear slightly misty after being sterilized with boiling water or vapor.

### Damaged accessories or buying of accessories

- Please go to the shop where the breast pump was purchased or find your local agency.

# Product Specification

## Content

Parts Name	Material	Heat Resistance	Place of Origin	Quantity
Electric Breast Pump	ABS Propylene Butadiene styrene copolymer	-20°C~80°C	Taiwan	1
Adaptor (100V~240V)				1
USB Power Cord				1
Silicone Diaphragm	Silicone	-20°C~120°C	Taiwan	2
Tube			Taiwan	1
Tube (T-Shaped)			Taiwan	1
Silicone Cover			Taiwan	2
Valve			Taiwan	2 + 2 (spare)
Storage Cap			China	2
Electric Breast Pump Top Cover	PP Polypropylene	-20°C~100°C	Taiwan	2
Breast Shield			Taiwan	2
Feeding Bottle Base			Taiwan	2

## Combi Standard Feeding Bottle

Parts Name	Material	Heat Resistance	Place of Origin	Quantity
Nipple	Silicone	-20°C~120°C	China	2
Feeding Bottle	PES Polyether Sulfone	-20°C~180°C	Taiwan	2
Nipple Ring	PP Polypropylene	-20°C~110°C	China	2
Hood			China	2

### Product Name :

Combi Natural Rhyme Double Electric Breast Pump

**Model Number :** DEBP17PI

**Net Weight :** Main Device 206 (±3) grams

### Dimensions :

**Packaging Dimension:** 38.3 x 26.5 x 11.2 cm

**Main Device :** 93.8 x 93.5 x 56.5 (±3) mm

**Rated Voltage (V) :** 100~240 VAC

**Rated Frequency (Hz) :** 50/60 Hz

**Total Rated electric Power Consumption :** 15W

**Manufacturing Date :** As printed on the packaging

**Manufacturing Code :** As printed on the base of product

**Place of Origin :** Taiwan

### Precautions or Warnings :

Please refer to the user manual for related information

### Instructions and Emergency Treatment :

Please refer to the user manual for related information